

Media for dechlorinators from the Löffler lab

<http://web.utk.edu/~microlab/LoefflerLab/Home.html>

Mineral Salts Medium ("DCB-1 medium")

Protocol

Salts	1 x [g/L]	100 x stock [g/L]
NaCl	1.0	100.0
MgCl ₂ x 6 H ₂ O	0.5	50.0
KH ₂ PO ₄	0.2	20.0
NH ₄ Cl	0.3	30.0
KCl	0.3	30.0
CaCl ₂ x 2 H ₂ O	0.015	1.5

To prepare medium	1 L
100 x salts	10 ml
Trace element solution	1 ml
Se/Wo solution	1 ml
Resazurin (0.1% solution)	0.25 ml
Electron donor	optional
H ₂ O bidest.	add 1000 ml

Some chemicals (e.g. fumarate, pyruvate) should not be autoclaved

- boil, cool down to room temperature under flushing with N₂/CO₂ (80/20)
- add 0.2 mM L-cysteine (MW 157.6) 0.031 g/L
- add 0.2 mM Na₂S x 9 H₂O 0.048 g/L [or: add 1 mM DL-dithiothreitol 0.154 g/L]
- add 30 mM NaHCO₃ 2.52 g/L
 adjust pH to 7.2 - 7.3 with CO₂ (to lower the pH with CO₂ is less time consuming)
- flush serum bottles with N₂/CO₂
- dispense medium (ca. 98 ml per bottle), close the bottles with black rubber stoppers (or Teflon-lined stoppers)
- autoclave after medium turned clear; autoclave bottles in a closed basket; do **NOT** remove bottles from autoclave before temperature of the medium has reached <70°C

Electron donors (formate, acetate, lactate, propionate) and acceptors (PCE or *cis*-DCE) can be added prior to heat sterilization. Fumarate, pyruvate are heat-labile and should not be autoclaved.

Trace element solution

Per liter: HCl (25% solution, w/w), 10 ml; FeCl₂ x 4 H₂O, 1.5 g; CoCl₂ x 6 H₂O, 0.19 g; MnCl₂ x 4 H₂O, 0.1 g; ZnCl₂, 70 mg; H₃BO₃, 6 mg; Na₂MoO₄ x 2 H₂O, 36 mg; NiCl₂ x 6 H₂O, 24 mg; CuCl₂ x 2 H₂O, 2 mg

Se/Wo solution

Per liter: 6 mg Na₂SeO₃ x 5 H₂O, 8 mg Na₂WO₄ x 2 H₂O and 0.5 g NaOH

Wolin Vitamins

Wolin, F. A., M. J. Wolin, and R. S. Wolfe. 1963. Formation of methane by bacterial extracts. J. Biol. Chem. **238**:2882-2886.

Vitamins	1000 x [mg/L]	Final conc.
biotin	20 mg/L	0.02 mg/L
folic acid	20 mg/L	0.02 mg/L
pyridoxine hydrochloride	100 mg/L	0.1 mg/L
riboflavin	50 mg/L	0.05 mg/L
thiamine	50 mg/L	0.05 mg/L
nicotinic acid	50 mg/L	0.05 mg/L
pantothenic acid	50 mg/L	0.05 mg/L
vitamin B12	1 mg/L	0.001 mg/L
p-aminobenzoic acid	50 mg/L	0.05 mg/L
thioctic acid	50 mg/L	0.05 mg/L

adjust pH to ~7.5 with 10 M NaOH (takes some time)

aliquot in 20 ml portions, freeze, store in dark place (light sensitive)

prepare a 200 x working stock solution, filtersterilize